

Replacement sterile filter cartridge It / min 400

Cod. TA-FT-001-C



Replacement sterile filter cartridge for TA-FT-001 It /min 400 to be replaced every 2 years.

The sterile MGF filters are placed at the air outlet of the compressor, the range of these filters covers all MGF dental compressors and have the following characteristics:

- The filter has a stainless steel cartridge that can be sterilized in an autoclave.
- It is recommended to sterilize the cartridge weekly for the best filtration efficiency.
- The filter can be installed with all oil-free compressors on the market.
- The filter is supplied with 2 drain taps to open it easily for removing the stainless steel cartridge.