



## MINI BOX 24/10 GENESI M

**Cod. D-MB150-024-M-15**

Oil free compressor, super silenced, FAD at 5 bar 120 l/min, tank capacity 24 L with membrane air dryer, filtration up to 0,01 micron and complete with MINI BOX acoustic cabinet. The MINI BOX acoustic cabinet provides an optimal noise reduction (50 dBA), the lowest noise level in this product category. User friendly and compact, it can be placed everywhere, including working environments.


 Maximum simultaneous use: 1

 Dental Chairs: 2

Medical Device CE 0425

										
Air Dryer	FAD <small>Fad @ 5bar</small>		Motor	Voltage Hz	Receiver	Noise Level	Duty Cycle	max P.	Size	net Weight
Yes	m³/min 0.135	l/min 135	kW 1.1	230V 50Hz	l 24	db(A) 50	50 %	bar 7	mm 530 x 500 x 800	Kg 67

## Air delivery performance



## Maintenance kits

Kit.	Description	Type	Hours
TA-FL-024-K-SCC	<p><b>GENESI M ordinary maintenance kit</b></p> <p>Ordinary maintenance kit for GENESI compressors with dryer. The kit consists of the cartridges of the dryer filtration battery (2 cartridges: 5 + 0.01 micron) and the relative O-rings, the pumping element suction filter and the non-return valve pad. Annual replacement of kit components is recommended. Maintenance operations can be performed at a higher frequency if the environment in which the compressor works is particularly dirty.</p>	Ordinary	1.000h/1y
MG-GC-011-K	<p><b>Extra ordinary maintenance kit for GENESI S and M versions</b></p> <p>The GENESI extraordinary maintenance spare parts kit includes: a suction filter, two elastic rings, O-rings, a line filter cartridge and two complete valve plates. The recommended maintenance intervals for this compressor are 1500 hours / 2 years. Maintenance operations can be performed at a higher frequency if the environment in which the compressor works is particularly dirty.</p>	Ordinary	1.500h/2y